

Colorimetry - To determine amount of nitrite present in the given water sample colorimetrically.



EXPERIMENT: 2

Date / /
Page 4

- Aim : To determine the amount of nitrite present in the given water sample colorimetrically.
- Requirements : Sodium nitrite (A.R.), sulphanilic acid, α -naphthyl amine, sodium acetate, HCl, pipettes, 100 cm³ standard measuring flasks, spectrophotometer, beakers, cuvettes, etc.
- Theory : The presence of nitrite ions (NO_2^-) in water sample may cause serious problems. It may react with dissolved oxygen affecting aquatic life. Therefore, the concentration of NO_2^- ion in water should be brought to a reasonable limit.

WHO suggested that the upper limit for it in drinking water is 80 ppm.

When a dilute acid solution of sulphanilic acid & α -naphthylamine are acted upon by NO_2^- ion or HNO_2 , a red colour is developed which has absorption maximum at 520 nm.



The sulphanilic acid is converted into corresponding diazo compound, which couples with α -naphthylamine to form a red coloured azo dye. So, the intensity of colour developed is directly proportional to the concentration of nitrogen (nitrite ion) in water, so it can be used to determine the amount of nitrite present in the water sample.

Teacher's Sign: _____

• Observations :

- 1) Prepare the reagent flask as given in the procedure.
- 2) Wavelength : 520 nm.

• Observation Table :

Flask No.	Volume of 10ppm of Nitrite sol ⁿ in cm ³	Sulphanilic acid sol ⁿ in cm ³	α -naphthyl amine sol ⁿ in cm ³	2M Sodium acetate sol ⁿ in cm ³	Conc. of nitrite ions in ppm	Absorbance at 5 nm A
1	0.0 (Blank)	2	2	2	-	-
2	1.0	2	2	2	0.1	0.13
3	2.0	2	2	2	0.2	0.22
4	3.0	2	2	2	0.3	0.27
5	4.0	2	2	2	0.4	0.29
6	5.0	2	2	2	0.5	0.30
7	1.0 (Sample)	2	2	2	"	A _s = 0.03

EXPERIMENT: No.

Date / / Page 5

• Procedure :

A] Preparation of standard solutions & reagents:

- 1) Standard nitrite solution: Dissolve 4.93 g of A.R. sodium nitrite in a minimum amount of distilled water. After shaking, make up to 1 dm³. This gives stock solution of nitrite ion. Dilute 10 cm³ of this solution to 1000 cm³ with distilled water in another standard measuring flask. This forms 10ppm of standard nitrite solution.
- 2) Sulphanilic acid reagent: Dissolve 0.60 g of A.R. sulphanilic acid in 100 cm³ of 20% hydrochloric acid in 100 cm³ standard measuring flask.
- 3) α -naphthyl amine reagent: Dissolve 0.48 g of A.R. α -naphthyl amine in 100 cm³ of 1.3% HCl (v/v) in 100 cm³ standard measuring flask.
- 4) 2M sodium acetate solution: Dissolve 16.4 g of A.R. anhydrous sodium acetate in about 60 cm³ distilled H₂O & make up the volume to 100 cm³ using 100 cm³ standard measuring flask.

B] Preparation of Calibration Curve:

- 1) Take 1, 2, 3, 4, 5 cm³ of 10ppm standard nitrite solution in five serially numbered 100 cm³ standard measuring flask.
- 2) Add 2 cm³ of sulphanilic acid solution to each flask. Shake well & wait for 5 minutes.
- 3) Now add 2 cm³ of α -naphthylamine followed by 2 cm³ of 2M sodium acetate to each flask. Wait for 10 minutes.
- 4) Dilute up to mark with fresh distilled water & shake.

Silabhi

Teacher's Sign



Shot on realme U1

• Calculations :

- 1) Absorbance of 100cm^3 of diluted water sample (A_s) = 0.03
- 2) Concentration of nitrite nitrogen of 100cm^3 of diluted water sample solution (x) = 0.04 ppm (From Graph 2.1)
i.e. $0.04\text{ mg} / 1000\text{cm}^3$


$\therefore 100\text{cm}^3$ of the diluted water sample solution contains 0.004 mg nitrite nitrogen.

i.e. 1cm^3 of water sample contains 0.004 mg nitrite nitrogen.

Now, 14 mg of Nitrogen = 46 mg of nitrite

$\therefore 0.004\text{ mg}$ of Nitrogen = $\frac{46 \times 0.004}{14}$

$$= \frac{0.184}{14} = 0.013\text{ mg of Nitrite (NO}_2\text{)}$$



• Calcul
 1) Abs
 2) Coni
 wat
 i.e
 ∴ 100
 0.1
 i.e. 10
 Now,
 ∴ 0.0

EXPERIMENT: No.

Date / /
Page 7

- 5) In another 100 cm³ standard measuring flask, prepare blank solution i.e. adding all reagent mentioned above without adding any amount of nitrite solution.
- 6) Switch on spectrophotometer & set the wavelength at 520 nm.
- 7) Measure absorbance (A) of each solution against blank prepared in step no 5.
- 8) Record the observations.

C] Analysis of sample solution :

- 9) Pipette out 1 cm³ of supplied water sample in another 100 cm³ standard measuring flask.
- 10) Add all the reagent in the same manner as in step no 2 & 3.
- 11) Dilute with distilled water & shake & wait for 5 minutes.
- 12) Determine absorbance of water sample against blank at 520 nm as reading.
- 13) Plot graph absorbance/O.D. against concentration of nitrite in ppm to form calibration curve.
- 14) With help of A & graph, determine the amount of nitrite nitrogen present in supplied water sample solution.

• Results:

- 1) Absorbance of 100 cm³ diluted water sample solution (A_s) = 0.03
- 2) Amount of nitrite in 1 cm³ of water sample = 0.0131 mg
- 3) Concentration of nitrite in the supplied water sample solution = 32.8×0.04 mg
= 13.12 mg

Teacher's Sign:

R

Shot on realme U1

Spectrophotometry: To determine amount of Fe(II) & Fe(III) in mixture using 1,10 phenanthroline spectrophotometrically.

EXPERIMENT: No. 3

Date / / Page 8

o Aim : Simultaneous determination of Fe(II) & Fe(III) by 1,10 phenanthroline method.

o Requirements :

- Standard Fe(II) solution: Dissolve 0.176 g of (AR) Ferrous ammonium sulphate in distilled water containing 2.5 cm³ of conc HCl. Dilute the solution with distilled water to 250 cm³. This gives standard Fe(II) solution of strength 100 ppm.
- 0.25% solution of 1,10 phenanthroline solution in distilled water.
- 10% Hydroxylamine hydrochloride solution.
- Dissolve 1.021 g of (AR) potassium hydrogen phthalate in 100 cm³ of distilled to prepare 0.05M potassium hydrogen phthalate solution.
- Sample solution containing a mixture of Fe(II) & Fe(III) in 100 cm³ standard measuring flask.
- 100 cm³ standard measuring flasks, distilled water, burette, 10 cm³ pipette, spectrophotometer, etc.

o Theory : Iron (II) reacts with 1,10 phenanthroline to form an orange-red coloured complex $[\text{C}_{12}\text{H}_8\text{N}_2]_3\text{Fe}^{2+}$. The colour intensity is independent of the acidity in the pH range 2-9 & is stable for long periods. It shows absorption at 515 nm at which iron (III) has no absorption.

Both Fe(II) & Fe(III) ions can be determined at 515 nm with pH of about 4 by reducing Fe(III) to Fe(II) which will give absorption due to both Fe(II) as well as Fe(III) .

R

Shot on realme U1

Sulaiman

Teacher's Sign:

• Observation :
Wavelength : 515 nm

• Observation Table :
Calibration Curve for Fe(II):

Flask No.	Standard Fe(II) Solution (cm ³)	Volume of 0.05M potassium hydrogen phthalate (cm ³)	Volume of 0.25% 1, 10 phenanthroline reagent (cm ³)	Concentration of Fe(II) ppm	Absorbance (A)
0	0.0	10	10	-	-
1	1.0	10	10	0.1	0.07
2	2.0	10	10	0.2	0.14
3	3.0	10	10	0.3	0.20
4	4.0	10	10	0.4	0.25
5	5.0	10	10	0.5	0.30
6	10cm ³ diluted sample solution	10	10	y	A = 0.12
7	10cm ³ diluted sample solution + 10cm ³ hydrogen amine solution	10	10	x	B = 0.28

• Observations :

From Graph 3.1,

x = Total amount of iron in [Fe(II) + Fe(III)] ppm

$x = 0.4$ ppm

y = Amount of Fe(II) in ppm

$y = 0.17$ ppm

$x - y$ = Amount of Fe(III) in ppm.

$0.4 - 0.17 = 0.23$ ppm

Shot on realme U1

EXPERIMENT: No.

Date / /
Page 9

While another aliquote without reduction of Fe(III) to Fe(II), will give absorption due to Fe(III) only. The difference between two gives amount of Fe(III) ion.

• Procedure :

I. Calibration Curve for Fe(II):

- Take eight 100cm³ standard measuring flasks serially numbered as 0, 1, ..., 7.
- Take 100ppm Fe(II) solution in a standard measuring flask No. 1 to 6 in quantity as given in Table.
- Add 10 cm³ of 0.05M potassium hydrogen phthalate to each flask.
- Add 10 cm³ of 0.25% 1, 10 phenanthroline reagent to each flask & dilute the solution up to mark with water.
- Mix the solution well & after 10 min, determine the absorbance at 515nm against reagent blank solution in flask No. 0.
- Plot the graphs of absorbance against concentration of solution at wavelength 515 nm.

II. Analysis of sample solution :

- Dilute the given sample of mixture of Fe(II) & Fe(III) to 100cm³ with distilled water in a standard measuring flask.
- To determine Fe(II) only:
 - Pipette out 10cm³ of diluted sample solution in 100cm³ standard measuring flask No. 7. Add 10cm³ 0.05M potassium hydrogen phthalate solution followed by 10cm³ of 0.25% 1, 10 phenanthroline solution.

Signature

Teacher's Sign

• Calculations :

From Graph 3.1,

1) Amount of Fe(II) present in 10cm^3 of diluted sample solution = 0.17 ppm

\therefore Amount of Fe(II) present in the given sample solution
 $= 10 \times 0.17 \text{ ppm} = 1.7 \text{ ppm}$ i.e. $1.7 \text{ mg}/1000 \text{ cm}^3$

$\therefore 100\text{cm}^3$ contains 0.17 ppm of Fe(II) ions.

2) Amount of Fe(III) present in 10cm^3 of diluted sample solution = 0.23 ppm.

\therefore Amount of Fe(III) present in the given sample solution
 $= 10 \times 0.23 \text{ ppm} = 2.3 \text{ mg}/1000 \text{ cm}^3$.

$\therefore 100\text{cm}^3$ contains 0.23 mg of Fe(III) ions.

EXPERIMENT: No.

Date / /
 Page 10

- ix) Dilute the solution to 100cm^3 with distilled water
- x) Measure the absorbance of solution after 10 min against reagent blank at 515nm . From the absorbance value, find out amount of Fe(II) from calibration curve (say $y \text{ ppm}$)
- b) To determine total Fe(II) & Fe(III) amount of iron:
 - xi) Pipette out 10cm^3 of diluted sample solution in 100cm^3 measuring flask No. 8
 - xii) Add 10cm^3 of hydroxylamine hydrochloride to it & shake
 - xiii) Then add 10cm^3 of 0.05M potassium hydrogen phthalate & 10cm^3 of 0.25% 1,10 phenanthroline reagent to it & shake well
 - xiv) Dilute the solution upto the mark with distilled water & shake well
 - xv) Measure the absorbance of solution at 515nm . From the absorbance value, find out total amount of Fe(II) & Fe(III) both from calibration curve (say $x \text{ ppm}$)
 - xvi) Therefore, amount of Fe(III) present can be calculated from $(x-y)$ values (say $z \text{ ppm}$).

Teacher's Sign

Graph : 3.1

Absorbance Vs. Conc. of Fe ion in ppm

Scale:


On X axis: 1cm = 0.05 ppm

On Y axis: 1cm = 0.0025

Amount of Iron in $[Fe(II) + Fe(III)]$ ppm = $x = 0.4$ ppm

Amount of $Fe(II)$ in ppm = $y = 0.17$ ppm

Amount of $Fe(III)$ in ppm = $z = x - y$
 $= 0.4 - 0.17$
 $z = 0.23$ ppm



o Results :

1) Amount of Fe(II) ion in a given mixture of Fe(II) & Fe(III) = $x = \underline{0.17 \text{ mg}}$.

2) Amount of Fe(III) ion in a given mixture of Fe(II) & Fe(III) = $z = \underline{0.23 \text{ mg}}$.

• Observation :

I] Standardisation of 0.1N HCl Solution :

Solution in burette = 0.1 N HCl

Solution in conical flask = 10 cm³ of diluted solution of borax

+ 15 cm³ water + 1-2 drop of 1% methyl orange indicator

Indicator = 1% methyl orange indicator

End point = Orange to pink.

• Observation Table :

Level	Burette Reading (cm ³)			C.B.R. (cm ³)
	I	II	III	
Initial	0.0	0.0	0.0	
Final	8.1	8.0	8.0	
Difference	8.1	8.0	8.0	$x = 8.0 \text{ cm}^3$

II] Estimation of Na₂CO₃ in Washing Soda :

Solution in burette = 0.1 N HCl

Solution in beaker = 10 cm³ of diluted solution of sample + 25 cm³ of distilled water.

Weight of the washing soda sample (W₁) = 0.5 g.

• Observation Table :

Obs. No.	Vol. of 0.1N HCl solution added in cm ³	Volume change (V ₂ - V ₁) in cm ³	Corrected Vol. (V _{cor}) = V + ΔV cm ³	pH of the solution	ΔpH	ΔV
1	0.0			10.40		
2	1.0	1	2	10.09	0.31	0.31
3	2.0	1	3	9.75	0.34	0.34
4	3.0	1	4	9.42	0.33	0.33
5	4.0	1	5	9.02	0.40	0.40
6	5.0	1	6	8.36	0.66	0.66
7	6.0	1	7	6.52	1.84	1.84

pH - metry : To determine percentage purity of sodium carbonate in washing soda pH

EXPERIMENT: 1


Date: / / Page: 1

• Aim : To determine percentage purity of sodium carbonate in washing soda, pH metrically.

• Requirements : Washing soda powder, 0.05N borax solution, 0.1N HCl solution, 1% methyl orange indicator, buffer solution of pH 4.01 & 9.18, burette, 100 cm³ standard measuring flasks, 100 cm³ beakers, etc.

• Theory : Na₂CO₃ is salt of a weak acid & a strong base. As a result, in aqueous solution, it hydrolyses to form a basic solution.

The strength & amount of Na₂CO₃ can be estimated with fair accuracy using pH metric titration technique against strong acid like HCl. The reaction is given as:



During the course of titration, pH is noted & graph of pH against volume of titrant added is plotted & from it, equivalence point is determined graphically. For accurate determination of equivalence point a first derivative graph of pH/V against the volume of titrant added is plotted.

• Procedure :


A] i] Preparation of sample solution :

- 1) Weigh accurately 0.5 g of given sample of washing soda & transfer it to a 100 cm³ beaker.
- 2) Dissolve in minimum quantity of distilled water. Then it to a 100 cm³ standard measuring flask & dilute it to the mark with distilled water.

Obs. No.	Vol. 0.1N HCl added in cm ³	$\Delta V = V_2 - V_1$	$V_{\text{cor}} = V + \Delta V$	pH	ΔpH	$\frac{\Delta \text{pH}}{\Delta V}$
8	7.0	1	8	5.84	0.68	0.68
9	8.0	1	9	5.40	0.44	0.44
10	9.0	1	10	5.04	0.36	0.36
11	10.0	1	11	4.63	0.41	0.41
12	11.0	1	12	3.50	1.13	1.13
13	12.0	1	13	1.50	2.0	2.0
14	13.0	1	14	1.16	0.34	0.34
15	14.0	1	15	0.99	0.17	0.17
16	15.0	1	16	0.85	0.14	0.14
17	15.5	0.5	16.0	0.80	0.05	0.1
18	16.0	0.5	16.5	0.75	0.05	0.1
19	16.5	0.5	17.0	0.72	0.03	0.06
20	17.0	0.5	17.5	0.67	0.05	0.1
21	17.5	0.5	18.0	0.64	0.03	0.06
22	18.0	0.5	18.5	0.62	0.02	0.04

• Calculations :

i) Normality of HCl solution :



$$N_1 V_1 = N_2 V_2$$

$$M_1 \times \frac{1}{1} = 0.05 \times 25$$

$$\therefore \text{Normality of HCl solution} = \frac{0.05 \times 25}{x}$$

$$A = \frac{0.05 \times 25}{8}$$

$$A = 0.15 \text{ N}$$

EXPERIMENT: No.

Date / /
Page 2

- A] i) Preparation of 0.05N borax solution :
- 1) Weigh accurately 0.381g of borax & transfer it to a clean dry 100cm³ beaker.
 - 2) Dissolve it in a minimum quantity of water, transfer it to a 100cm³ standard measuring flask & dilute it upto mark with distilled water.
- B] Standardisation of 0.1N HCl solution :
- 3) Rinse & fill a burette with the supplied 0.1N HCl.
 - 4) Pipette out 0.5cm³ of 0.05N borax in clean 150cm³ conical flask. Add 2 drops of 1% methyl orange & about 10-15 cm³ distilled water to it.
 - 5) Titrate it against 0.1N HCl from burette. End point will be from yellow to pink colour (x cm³).
 - 6) Repeat & obtain two more constant readings.
- C] Estimation of Na₂CO₃ in washing soda sample :
- 7) Pipette out 10cm³ of washing soda sample solution in 100cm³ beaker. Add about 25 cm³ of distilled water to it.
 - 8) Standardise pH meter at pH 9.18 & 4.01 using buffer solution of borax & potassium hydrogen phthalate resp.
 - 9) Wash & dry electrodes & dip in sample solution & record the initial pH of the solution.
 - 10) Add 1 cm³ of 0.1N HCl from burette, stir & record the pH.
 - 11) Continue the addition of HCl till two distinct sharp pH are observed.
 - 12) Add HCl dropwise & till pH reaches to value about 2.
 - 13) Plot graphs of (i) pH Vs. volume of HCl added in cm³ (ii) $\Delta \text{pH} / \Delta V$ Vs. mean volume of titrant added.

Explain

Teacher's Sign

R

Shot on realme U1

II] Amount of Na_2CO_3 in the given washing soda sample :
From Graph,

i) pH against volume of HCl solution added :

$$V_1 = 5.6 \text{ cm}^3 \quad V_2 = 11.2 \text{ cm}^3$$

ii) ΔpH against mean volume of titrant added :

$$V_1 = 6 \text{ cm}^3 \quad V_2 = 12 \text{ cm}^3$$

\therefore Volume of HCl solution required for complete titration
= $V_2 = 11.6 \text{ cm}^3$

$$\begin{aligned} 2 \text{ moles of HCl are required to neutralise 1 mole of NaHCO}_3 \\ \therefore 2000 \text{ cm}^3 \text{ of } 1 \text{ N HCl solution} &= 106 \text{ g of Na}_2\text{CO}_3 \\ \therefore \text{cm}^3 \text{ of } 0.1 \text{ N HCl solution} &= \frac{106 \times 11.6 \times A}{2000} \\ &= \frac{53 \times 11.6 \times 0.15}{1000} \\ &= 0.092 \text{ g of Na}_2\text{CO}_3 \end{aligned}$$


10 cm³ of the diluted sample solution contains 0.092 g of Na_2CO_3 .

\therefore 100 cm³ of the diluted sample solution contains $10 \times 0.092 \text{ g}$ of $\text{Na}_2\text{CO}_3 = 0.92 \text{ g}$ of $\text{Na}_2\text{CO}_3 = Z \text{ g}$

Percentage of Na_2CO_3 in washing soda sample :

$$\begin{aligned} \% \text{ Na}_2\text{CO}_3 &= \frac{\text{Weight of Na}_2\text{CO}_3 \text{ observed}}{\text{Weight of washing soda sample}} \times 100 \\ &= \frac{0.092}{0.5} \times 100 = 0.184 \times 100 = 18.4\% \end{aligned}$$


Shot on realme U1



Graph 1

On X axis 5.1 cm = 1 cm³
 On Y axis 5.1 cm = 0.2 cm³

DP
 DV



- 14) Determine the volumes of V_1 & V_2 from graph corresponding to half & complete neutralisation.
- 15) From the volume V_2 , calculate the amount of Na_2CO_3 in washing soda sample.

o Results :

- 1) 10cm^3 of diluted sample solution required $(x) = 80$ of 0.15N HCl solution
- 2) Percentage of Na_2CO_3 present in washing soda sample = 18.4%

~~Supong~~
~~06/06/23~~

Observation Table:

Conc. of Cr (M)	Volume of 0.01 M $K_2Cr_2O_7$	Absorbance at 440 nm	Absorbance at 525 nm
0.01	10	0.416	0.055
0.0005	5	0.200	0.017
	Unknown	0.323	0.736

Conc. of Mn (M)	Volume of 0.01 M $KMnO_4$	Absorbance at 440 nm	Absorbance at 525 nm
0.00055	5.0	0.025	0.923
0.00025	2.5	0.014	0.277

Calculations:

I] Preparation of 0.1 M $K_2Cr_2O_7$

$$\text{Weight} = \frac{\text{mol. wt.} \times \text{conc} \times \text{volume}}{1000}$$

$$= \frac{294.19 \times 0.01 \times 250}{1000}$$

$$= 0.735 \text{ gm}$$

II] Volume for 0.01 M $K_2Cr_2O_7$

i) $M_1V_1 = M_2V_2$

$$0.01 \times V_1 = 0.001 \times 100$$

$$V_1 = \frac{0.001 \times 100}{0.01}$$

$$V_1 = 10 \text{ cm}^3$$

ii) $M_1V_1 = M_2V_2$

$$0.01 \times V_1 = 0.0005 \times 100$$

$$V_1 = \frac{0.0005 \times 100}{0.01}$$

$$V_1 = 5 \text{ cm}^3$$

Determine amount of Cr (VI) & Mn (VII) in given solution by simultaneous spectrophotometrically.

EXPERIMENT No. 4 Date / / Page 13

Aim: Simultaneous determination of Chromium & Manganese in a mixture by spectrophotometrically.

Requirements:

- 0.005 M $K_2Cr_2O_7$ standard solution: Prepare 0.1 M $K_2Cr_2O_7$ solution by dissolving 2.942 g of $K_2Cr_2O_7$ in 100 cm³ standard measuring flask containing 10 cm³ of 10 M H_2SO_4 & dilute to 100 cm³ distilled water. Stock solution $K_2Cr_2O_7$.
- Dilute 5 cm³ stock solution to 100 cm³ with 1 M H_2SO_4 in standard measuring flask. It gives standard 0.005 M $K_2Cr_2O_7$.
- 0.005 M $KMnO_4$ solution: Dissolve 0.315 g of $KMnO_4$ in 100 cm³ standard measuring flask containing 10 cm³ of 10 M H_2SO_4 & dilute to 100 cm³ with distilled water. This gives the stock solution of $KMnO_4$. Dilute 5 cm³ of stock $KMnO_4$ with 1 M H_2SO_4 upto 100 cm³. This gives 0.005 M standard $KMnO_4$.
- 1 M H_2SO_4 solution
- Sample solution containing Chromium (VI) & Manganese (VII) in 100 cm³ standard measuring flask.
- 6 100 cm³ standard measuring flasks, burette, 5 cm³ & 10 cm³ pipette, spectrophotometer, etc.

Theory:

When there is no interaction between two solute in a solution it is possible to determine their conc. in mixture.

R

Shot on realme U1

Sunam

Teacher's Sign

④ For Mn⁵²⁵

$$1) E_{Mn}^{525} = \frac{0.943}{0.0005} = 1886$$

$$ii) E_{Mn}^{525} = \frac{0.277}{0.00025} = 1108$$

$$\text{Average } E_{Mn}^{525} = \frac{1886 + 1108}{2} = 1497 \text{ cm}^{-1} \text{ mol}^{-1} \text{ dm}^{-3}$$

$$\therefore C_{Cr} = \frac{A_{mix}^{440} \cdot E_{Mn}^{525} - A_{mix}^{525} \cdot E_{Mn}^{440}}{E_{Cr}^{440} \cdot E_{Mn}^{525} - E_{Cr}^{525} \cdot E_{Mn}^{440}}$$

$$= \frac{(0.323 \times 1497) - (0.736 \times 53)}{(408 \times 1497) - (44.5 \times 53)}$$

$$C_{Cr} = 7.30 \times 10^{-4} \text{ mol/dm}^3$$

$$C_{Mn} = \frac{A_{mix}^{525} \cdot E_{Cr}^{440} - A_{mix}^{440} \cdot E_{Cr}^{525}}{E_{Mn}^{525} \cdot E_{Cr}^{440} - E_{Mn}^{440} \cdot E_{Cr}^{525}}$$

$$= \frac{(0.736 \times 408) - (0.323 \times 44.5)}{(1497 \times 408) - (53 \times 44.5)}$$

$$C_{Mn} = 4.70 \times 10^{-4} \text{ mol/dm}^3$$

For Mn, 1M KMnO₄ = 54.94 g of Mn

$$C_{Mn} = 4.70 \times 10^{-4} \text{ M KMnO}_4$$

$$C_{Mn} = (54.94) \times (4.70 \times 10^{-4})$$

$$= 25.82 \times 10^{-3} \text{ g of Mn}$$

$$= \underline{25.82 \text{ mg of Mn}}$$

EXPERIMENT: No.

Date / /
Page 15

i) Measure the absorbance at 420 nm & 520 nm using 1M H₂SO₄ as blank.

ii) Determination of conc. of chromium (Cr) & manganese (Mn) in mixture.

i) Dilute sample upto the mark with 1M H₂SO₄ in 100 cm³ standard measuring flask.

ii) Transfer 10 cm³ of diluted solution to another 100 cm³ standard measuring flask No. 3 & dilute upto the mark with 1M H₂SO₄.

iii) Measure absorbance at 420 nm & 520 nm using 1M H₂SO₄ as blank.

Sunsh

Teacher's Sign



Shot on realme U1

For Cr, $1M\ K_2Cr_2O_7 = 2 \times \text{At. wt of Cr}$
 $= 2 \times 52$
 $= 104$

$$\therefore C_{Cr} = 7.30 \times 10^{-4} M\ K_2Cr_2O_7 = 104 \times 7.30 \times 10^{-4}$$
$$= 0.07592\ g\ of\ Cr$$
$$= \underline{75.92\ mg\ of\ Cr}$$

EXPERIMENT:

No.

Date / /
Page 16

• Results:

- i) Amount of Cr(VI) present in the given solution = 75.92 mg
- ii) Amount of Mn(VII) present in the given solution = 25.82 mg

~~Supervisor~~
06/06/23

Dalal

Teacher's Sign.



Shot on realme U1

Volume for 0.01 M KMnO_4

$$i) M_1 V_1 = M_2 V_2$$

$$0.01 \times V_1 = \frac{0.0005 \times 100}{0.01}$$

$$V_1 = 5 \text{ cm}^3$$

$$ii) M_1 V_1 = M_2 V_2$$

$$0.01 \times V_1 = 0.00025 \times 100$$

$$V_1 = \frac{0.00025 \times 100}{0.01}$$

$$V_1 = 2.5 \text{ cm}^3$$

ii] $E = A/Ic$ where A is absorbance

C = concentration

I = path length

The molar extraction coefficient as $C_{Cr^{4+0}}$, $M_{Mn^{4+0}}$, $C_{Cr^{5+25}}$, $M_{Mn^{5+25}}$ calculate by

$$\lambda_{max}^{440} = E_{Cr^{4+0}} C_{Cr^{4+0}} + E_{Mn^{4+0}} C_{Mn^{4+0}} \quad \text{--- (1)}$$

$$\lambda_{max}^{525} = E_{Cr^{5+25}} C_{Cr^{5+25}} + E_{Mn^{5+25}} C_{Mn^{5+25}} \quad \text{--- (2)}$$

i) For Cr^{4+0}

$$i) E_{Cr^{4+0}} = \frac{0.416}{0.001} = 416$$

$$ii) E_{Cr^{4+0}} = \frac{0.2}{0.0005} = 400$$

$$\text{Average } E_{Cr^{4+0}} = \frac{416 + 400}{2} = 408 \text{ cm}^{-1} \text{ mol}^{-1} \text{ dm}^{-3}$$

ii) For Cr^{5+25}

$$i) E_{Cr^{5+25}} = \frac{0.065}{0.001} = 65$$

$$ii) E_{Cr^{5+25}} = \frac{0.017}{0.0005} = 34$$

$$\text{Average } E_{Cr^{5+25}} = \frac{65 + 34}{2} = 49.5 \text{ cm}^{-1} \text{ mol}^{-1} \text{ dm}^{-3}$$

iii) For Mn^{4+0}

$$i) E_{Mn^{4+0}} = \frac{0.025}{0.0005} = 50$$

$$ii) E_{Mn^{4+0}} = \frac{0.014}{0.00025} = 56$$

$$\text{Average } E_{Mn^{4+0}} = \frac{50 + 56}{2} = 53 \text{ cm}^{-1} \text{ mol}^{-1} \text{ dm}^{-3}$$

The total absorbance of solution at given wavelength is equal to sum of absorbance of the individual components.

$$A_{mix}^{\lambda_1} = A_1^{\lambda_1} + A_2^{\lambda_1} \quad \text{--- (1)}$$

$$A_{mix}^{\lambda_2} = A_1^{\lambda_2} + A_2^{\lambda_2} \quad \text{--- (2)}$$

Absorbance is directly proportional to the product of molar absorptivity (E) & conc C . Hence for path length of 1 cm, we have

$$A_{mix}^{\lambda_1} = E_1^{\lambda_1} C_1 + E_2^{\lambda_1} C_2$$

$$A_{mix}^{\lambda_2} = E_1^{\lambda_2} C_1 + E_2^{\lambda_2} C_3$$

Solution of the simultaneous equation gives

$$C_1 = \frac{(E_2^{\lambda_1} A_{mix}^{\lambda_2}) - (E_2^{\lambda_2} A_{mix}^{\lambda_1})}{(E_1^{\lambda_1} A_{mix}^{\lambda_2}) - (E_1^{\lambda_2} A_{mix}^{\lambda_1})}$$

$$C_2 = \frac{(E_1^{\lambda_1} A_{mix}^{\lambda_2}) - (E_1^{\lambda_2} A_{mix}^{\lambda_1})}{(E_2^{\lambda_1} A_{mix}^{\lambda_2}) - (E_2^{\lambda_2} A_{mix}^{\lambda_1})}$$

Procedure:

i) Determination of molar absorptivity at 440 nm & 525 nm

ii) Take four numbered 100 cm³ standard measuring flasks

3) In flask No. 1 & 2 add 10 cm³ & 20 cm³ standard potassium dichromate solution & dilute upto the mark with 1M H_2SO_4 . This gives 0.0005 M & 0.001 M potassium dichromate solution.

4) In flask No. 3 & 4 add 5 cm³ of standard $\text{K}_2\text{Cr}_2\text{O}_7$ dilute upto mark with 1M H_2SO_4 . This gives 0.00025 M & 0.0005 M $\text{K}_2\text{Cr}_2\text{O}_7$.